#### INSTRUCTIONS

Thank you for purchasing the Vision Assessment Corporation Random Dot 3s Stereo Acuity Test, P/N 1009s.



#### **PURPOSE**

This test uses a random dot pattern for testing extra fine stereo discrimination (depth perception) to rapidly test for amblyopia, strabismus and other binocular acuity problems in adults, early & non-readers and non-verbal children & adults.

Depth Perception: (n) The ability to perceive spatial relationships, especially distances between objects, and three dimensions.


Geometric Shapes: Triangle, Circle, Square & Plus

Δ	Triangle
0	Circle
	Square
C)	Plus

## **FAMILIARIZE YOURSELF WITH THE TEST**

- Test consists of:
  - 1. 1 Booklet 2 Vectographic Pages
    - Page 1

4 levels of gross disparity




Page 2 (Section A)
10 levels of disparity for critical testing



Page 2 (Section B)
3 gross levels of disparity



2. 1 pair adult polarized viewers



3. 1 pair pediatric polarized viewers



4. Answer Key on back cover



#### **TESTING CONDITIONS**

- Well-lit, glare-free area
- If reflections or glare on the vectographic pages can be seen, try tilting the booklet or choose another testing location

# **ADMINISTERING THE TEST**

- Page 1
  - 1. Place the polarized viewers on the test subject
  - 2. **PLEASE NOTE**: IF TEST SUBJECT WEARS PRESCRIPTION GLASSES, THE POLARIZED VIEWERS SHOULD BE WORN OVER GLASSES
  - 3. Before opening test booklet, familiarize test subject with the 4 Shapes on front cover (*Triangle, Circle, Square & Plus*)
  - 4. The test subject should be able to correctly identify at least 3 of the 4 non-stereo, two-dimensional Random Dot shapes on the front cover before administering the test
  - Open test booklet and fold so that only Page 1 can be viewed by test subject
  - 6. Hold the test approximately 16 inches (40 cm) from the test subject
  - 7. Each Box A-D contains a floating Shape (*Triangle, Circle, Square or Plus*), each with decreasing levels of Stereopsis. The Shape in Box A also acts as floating Demonstration Target to help familiarize the test subject with the test
  - 8. Starting with Box A, ask the test subject to identify the Shape
  - 9. Repeat previous step for boxes B, C & D
  - 10. If test subject is unable to identify at least 1 Shape or is only able to correctly identify the Demonstration Shape in box A, turn to Page 2 and ask the test subject to identify the "L" & "R" in the center of the page. If the patient is unable to identify both the "L" & "R", the patient is untestable

## Scoring

- 1. Refer to Answer Key on back of test booklet
- 2. Record the level of stereopsis for the last Shape correctly identified
- Page 2, Section A (For Adults)
  - 1. Turn booklet so that only page 2 can be viewed by test subject
  - 2. Ask test subject to look at the 10 Circle targets each containing 4 smaller Circles

- 3. Only one Circle in each target will appear to be floating forward, and the other 3 Circles will appear to be recessed or floating backwards
- 4. Starting with Circle #1, ask the test subject to point to or identify (top, bottom, left or right) the Circle that appears to be floating forward or different from the others
- 5. Repeat last step for Circles 2-10
- 6. When one is target is identified incorrectly, retest the preceding target to ensure that the test subject can achieve this level and was not just guessing

## Scoring

- 1. Refer to Answer Key on back of test booklet
- 2. Record the level of stereopsis for the last target in which the test subject correctly identified the floating Circle

## Section B (For Children)

- 1. Each Row A, B & C consists of 4 boxes each containing a Shape that floats forward
- 2. The Shapes in each Row are set at the same level of stereopsis with Row A at 400 seconds, Row B at 200 seconds & Row C at 100 seconds of arc
- 3. Starting with Row A, ask the test subject to identify the Shape in each box (Triangle, Circle, Square or Plus)
- 4. Repeat previous step for Rows B & C

#### Scoring

- 1. Refer to Answer Key on back of test booklet
- 2. The test subject must correctly identify 3 of the 4 Shapes in a given Row to achieve that level of stereopsis
- 3. Record the level of stereopsis of the last Row for which the test subject correctly identified at least 3 of the 4 Shapes


SCORE KEY	Seconds of arc 16 in., (40 cm)	Shepard Percentage	Verhoff Distance
Α	400	15%	0.1
В	200	30%	0.2
С	100	50%	0.3

# **CARE/HANDLING & STORAGE**

- Clean vectographic pages with a soft, damp lint-free cloth
- \* Store closed booklet in a cool, dry place away from direct sunlight
- Clean polarized viewers using lens cleaner and soft, lint free cloth

#### WARRANTY

• 1 year manufacturer warranty from date of purchase



# Vision **Assessment** Corporation

Random Dot 3 s Stereo Acuity Test

P/N 1009s **INSTRUCTIONS** 

**Vision Assessment Corporation** 2675 Coyle Avenue Elk Grove Village, Illinois 60007 USA

Phone: 1 847 239 5889

Toll Free USA: 1 866 887 9692

Email: sales@VisionAssessment.com Web: www.VisionAssessment.com

