

Optimum Convenience for Maximum Competence:

OT-2001

OPHTHALMIC
OPERATION CHAIRS




OT-2001

OPHTHALMIC OPERATION CHAIRS

**Maximum-Comfort
Mobile Recliner Chairs
For Operations, Treatment & Transportation.**

TECHNO-SPECIALITIES:

- Custom-developed to satisfy ALL requirements of ophthalmologists.
- The patient can be efficiently prepared for the operation, operated on, and offered post-operative care on, OT-2001.
- Easy transformation from a smooth transportation chair into a versatile operation table.
- Ergonomic design, with various position settings.
- Multi-Adjustable Head Rest can be tilted independently and smoothly height-adjusted, enabling the patient's head to be kept in a convenient position.
- Ample leg-room for surgeons.
- High load-bearing, stable Twin Wheel Castors with direction-locking.
- Hand control panel and foot control for all movements.
- Smooth, CE-marked German DC Motor-controlled chair lifting and back & forth positioning.



Specifications	
Height range min. / max.	580 mm / 810 mm
Total length	1860 mm
Total width	610 mm
Lifting speed	8 mm/s
Backrest section	Electro-cylinder
Weight capacity	135 Kgs.
Input tension	230 V 50/60 Hz. 110 V 50/60 Hz.
Gross weight	118 Kgs.