

Neitz Binocular Loupes BL series

Headband Type



NEITZ

Made by Neitz for Your Medical Field

NEITZ Headband Loupes

Excellent balance and comfort for all medical professionals.

- Hinge joints support fast & easy vertical adjustment.

- Multiple adjustment points for desirable fit.

- Flip-up Function: easy to see unmagnified view



- **Adjustable PD & working distance**
Revolvable barrels enable you to adjust the working distance for your perfect vision. User-oriented unique function, made by NEITZ.

- **LED Spot illuminators** with outstanding brightness and illumination range are available. (optional)



NSI-III



NSI-X



Yellow Filter
(Optional for NSI-III)

Specification

| Model | Magnification | Working Distance | Field of View Diameter | Weight (incl. Headband) |
|-------------------------------------|---------------|--|--------------------------------------|-------------------------|
| Galilean Loupes | | | | |
| BLS-1 HB | 1.5x – 2.0x | 550 mm – 350 mm (21.7 in – 13.8 in) | 145 mm – 120 mm (5.7 in – 4.7 in) | 251 g |
| BLS-2 HB | 2.0x – 2.5x | 550 mm – 350 mm (21.7 in – 13.8 in) | 113 mm – 58 mm (4.5 in – 2.3 in) | 246 g |
| BLS-3 HB | 2.5x – 3.0x | 550 mm – 350 mm (21.7 in – 13.8 in) | 78 mm – 38 mm (3.1 in – 1.5 in) | 240 g |
| BLD-3 HB | 2.5x – 3.0x | 350 mm – 250 mm (13.8 in – 9.8 in) | 80 mm – 50 mm (3.2 in – 2.0 in) | 235 g |
| Prismatic (Keplerian) Loupes | | | | |
| BLP-4 HB | 4.0x – 4.5x | 550 mm – 350 mm (21.7 in – 13.8 in) | 80 mm – 50 mm (3.2 in – 2.0 in) | 258 g |
| BLP-6 HB | 5.5x – 6.0x | 550 mm – 350 mm (21.7 in – 13.8 in) | 40 mm – 25 mm (1.6 in – 1.0 in) | 261 g |