

Products details







Warning

High temperature alert

When the temperature is higher than 37.5°C , a "HI" will appear on the display screen of K3 infrared thermometer and a beep will be emitted to prompt timely attention to personal health.



Accuracy

High accuracy

K3 infrared thermometer can measure human body temperature in short time, the error accuracy is controlled between 0.1 and 0.2°C



Product introduction



Production Specification

Temp.measurement	Error +_0.2℃	Reaction time	0.5s
Type of alarm	High temp auto alert,shining and buzzing	transmission interface	USB standard interface
Measure distance	5-10cm auto measurement	Input voltage	USB DC 5V
Stand-by time	about one week	Product Net.Weight	350g
Installation	Tripod fixed,sponge tape fixed,wall hanging hole fixed.	Product Meaurement	170*115*140mm
Working temperature	10-40℃ (recommendate 15-36 ℃)	Measurement range	0-50 ℃
Standby mode	5s standby,red dot dynamic cycle	Display screen	Digital dispaly

Energy saving

Energy saving mode

Super long standby time, energy saving and environmental protection. Charge one time and use for a long time.



K3

Accurate temperature measurement,
one second reaction time

Easy operation




Multiple fixed mode

3 kinds of fixed mode fit for different using enviropment.



Tripod fixed



Sponge tape fixed



Wall nail hanging hole fixed

Smart infrared thermometer

SPEED

0.5s reaction time

Non-contact infrared thermometer,easy and fast operation,avoid cross infection.



PURPOSE

K3 infrared thermometer can set in crowded place,
just like:

railway station, bus station, hospital, company, school etc.



Packing details& information



Packing list

Product G.W	500g
Product N.W	350g
Product Size	170×115×140mm
Gift Box Size	205x150x177mm
Outer Carton size	460x370x425mm
Package per Carton	12 PCS