

SLIT LAMP (2 STEPS)

O-SL-01

Objective	: 1x, 6x
Total Magnification	: 16x, 25x
Slit Rotation	: 0-180°
Eye Piece Magnification	: 10x, 16x
Illumination Type	: Halogen
Slit Aperture	: 0.2mm, 2mm, 6mm, 10mm, 12mm
Color	: White
Halogen Lamp	: 06 Volts, 20 Watt
Slit Diaphragm Disc	: Wedge Shaped Diaphragm of Infinitely Variable Slit Lengths
Filters	: Cobalt Filter, Green Filter, Yellow Filter, Natural Density & Open Aperture

HIGH QUALITY
MEDICAL EQUIPMENTS BY

