



American  
Association for

Pediatric Ophthalmology  
and Strabismus

## Vision Screening Kits



### AAPOS Screening Kit (LEA SYMBOLS®)

The AAPOS Vision Screening Kit is for screening the vision of individuals beginning at age 3 through adulthood. Spiral-bound, multisided, distance visual acuity eye charts that can be hung on a wall or held in the hands. The distance charts contains both Full-Chart Threshold and Single Critical Line formats with the following optotypes:

- **Sloan Letters** for when children know their letters.
- **LEA SYMBOLS®** for preschoolers and for children who do not yet know their letters.

The AAPOS Vision Screening Kit includes Black Handheld Flip Paddle Occluder, Children's Flip Up Occluder Glasses, 10-foot cord for measuring, Translucent Response Key, Instructions and DVD with additional information.

120



### AAPOS Screening Kit (HOTV)

The AAPOS Vision Screening Kit is for screening the vision of individuals beginning at age 3 through adulthood. Spiral-bound, multisided, distance visual acuity eye charts that can be hung on a wall or held in the hands. The distance charts contains both Full-Chart Threshold and Single Critical Line formats with the following optotypes:

- **Sloan Letters** for when children know their letters.
- **HOTV** for preschoolers and for children who do not yet know their letters.

The AAPOS Vision Screening Kit includes Black Handheld Flip Paddle Occluder, Children's Flip Up Occluder Glasses, 10-foot cord for measuring, Response Key, Instructions and DVD with additional information.

129



### AAPOS Supplemental Screening Kit

The AAPOS Vision Screening Kit: Supplemental Screening Package is for screening the vision of individuals beginning at age 3 through adulthood.

The near vision chart contains both Full-Chart Threshold and Single Critical Line formats with the following optotypes:

- **Sloan Letters** for when children know their letters.
- **LEA SYMBOLS®** for preschoolers and for children who do not yet know their letters.

The AAPOS Supplemental Vision Screening Kit includes Black Handheld Flip Paddle Occluder, Flip Up Occluder Glasses, 16 Plate Color Vision Test, Random Dot Butterfly stereopsis test, Instructions and DVD with additional information.

130

### AAPOS Near Vision Kit

The AAPOS Vision Screening Kit is for screening the vision of individuals beginning at age 3 through adulthood. Spiral-bound, multisided, near eye charts with attached 16-inch cord. The near cards contains both Full-Chart Threshold and Single Critical Line formats with the following optotypes:

- **Sloan Letters** for when children know their letters.
- **LEA SYMBOLS®** for preschoolers and for children who do not yet know their letters.

The AAPOS Vision Screening Kit includes Black Handheld Flip Paddle Occluder, Children's Flip Up Occluder Glasses, 10-foot cord for measuring, Translucent Response Key, Instructions and DVD with additional information.

136



### AAPOS Near/Distance Vision Screening Kit

Kit Includes: School Nurse Vision Screening Tutorial DVD created by Kathy Lee, MD, PhD, a pediatric ophthalmologist

4 Double-Sided, 10-Foot Distance Eye Charts and 4 Double-Sided, 16-Inch Near Vision Charts with

- **LEA SYMBOLS®**
- **LEA NUMBERS®**
- **HOTV**
- **Sloan Letters**

Double-sided Response Panel with LEA SYMBOLS® and HOTV for Matching, Frosted Flip-Up Occluder Glasses.

110