

REFRACTION SYSTEM
CRT-8800



Power Arm and user-friendly Keypad allow you to experience more advanced refraction environment.

Power Arm

A smooth operating power arm with multi control sequence enhances the operator's convenience and efficiency.
(Left/Right Movement Angle : 45°)



Refraction chair

Unit chair, built on the ergonomic design, improves the satisfaction to the examinee, and the maximum convenience to the operator.

Keypad

Highly refined Key Pad with one-touch concept provides more convenient refraction environment at your fingertips.

Chair Stop Bar

Sensor-Operated Chair Stop Bar prevents the interference between the patient's legs and the table, allowing the comfortable space for the patient.

SPECIFICATION

Dimensions / Weight	1,035(W) x 970(D) x 1,652(H)mm / 140kg
Refractor Arm	Automated Power Arm Stroke : 400mm
Trial Lens Drawer	566(W) x 380(D) x 65(H)mm
Table Size	1,010(W) x 470(D) x 30(T)mm
Table Rotating Angle	90°
Head Lamp	50W
Power Supply	AC 110/220V, 50/60Hz
Power Consumption	420W (Max)
Others	Chair included (up/down stroke :150mm)

Designs and details can be changed without prior notice for improvements.