

FEATURES

- Simple but multifunction, with phoropter arm and projector stand
- Phoropter arm can be moved up and down, and turn around by manual
- Sliding table makes it easy to move the refractometer.

SPECIFICATIONS

Unit size.....	900x520x1710mm
Table size.....	900x500mm
Sliding plate size.....	300x400mm
Drawer size.....	543x370x50mm
Phoropter arm up and down range.....	190mm
Height of projector stand.....	1510mm
Power of lamp.....	110V/220V, 60HZ/50HZ, 15W
Power of unit.....	110V/220V, 60HZ/50HZ, 80W
Packing / Weight.....	980x590x420mm/60kgs

