

## CP-180A COMBINATION TABLE

WITH SLIDING PLATE

### FEATURES

- Simple but multifunction, with phoropter arm and projector stand
- Phoropter arm can be moved up and down, and turn around by manual
- Sliding table makes it easy to move the refractometer.

### SPECIFICATIONS

|                                      |                           |
|--------------------------------------|---------------------------|
| Unit size.....                       | 900x520x1710mm            |
| Table size.....                      | 900x500mm                 |
| Sliding plate size.....              | 300x400mm                 |
| Drawer size.....                     | 543x370x50mm              |
| Phoropter arm up and down range..... | 190mm                     |
| Height of projector stand.....       | 1510mm                    |
| Power of lamp.....                   | 110V/220V, 60HZ/50HZ, 15W |
| Power of unit.....                   | 110V/220V, 60HZ/50HZ, 80W |
| Packing / Weight.....                | 980x590x420mm/60kgs       |



CE FDA

## CP-180AT COMBINATION TABLE AND CHAIR

WITH SLIDING PLATE

Full functions of CP-180A, with electric chair  
Chair up and down range: 150mm



CE FDA

## CP-180B COMBINATION TABLE

WITH SLIDING PLATE

Full functions of CP-180A, but without projector stand



CE FDA