



ZN-185
ZN-215

LED VISUAL CHART
LED VISUAL CHART

18.5"
21.5"HD



ZN-185

CE FDA



ZN-215

CE FDA

Both table stand and wall mount are as standard accessories

SPECIFICATIONS

LED monitor.....ZN-185:18.5-inch TFT Panel, 1366x768 pixels
ZN-215:21.5-inch TFT Panel, 1920x1080 pixels
Brightness.....250CD/m²
Optotype.....E, C(Landolt),Alphabet, Numbers, Children
Visual unit.....Decimal,Snellen(ft,20/20),Snellen(m,6/6)
Distance.....1M~6M (in steps of 0.1m)
Mask type.....Single mask, Horizontal mask, Vertical mask,
Red/Green mask, White/Black mask,
Contrast adjustment, Red/Green adjustment
Special test masks.....Ishihara, Cross cylinder, Astigmatism,
Cross Grid, Fixation point Test, Red/Green Panel,
Worth Test, Red/Green Balance, Traffic light,
Coincidence Test, Stereo Test, Schober Test, ETDRS
Other functions.....Randomization Optotypes Playing function,
R/G setting function, Video playing function,
Screen-saver function, standby function

Standard accessories



Wall mount



Power adapter



Remote
