

SPECIFICATIONS

Magnification of Eyepiece.....	15X
Total Magnification.....	18X
Illumination.....	2X 6V/12W halogen bulb
Measurement Range:	
Radius of Curvature.....	5.5~11mm
Diopter.....	30~60Dpt
Reading Increment(radius).....	0.02mm
Measuring Accuracy.....	0.25Dpt
Compensation Refractive index n.....	1.3375
Minimum Surface required.....	r=5.5 Φ 1.65 r=7.5 Φ 2.3
For measurement.....	r=11 Φ 4.00(mm)
Angle Measuring Disc.....	0°~180°
Input Voltage.....	110V/220V

