

D220

LED VISUAL CHART

23" HD

POLARIZED SCREEN ECONOMICAL

- Big size monitor, clear display
23" TFT-LED monitor, 1920x1080 pixels, 250CD/m²
- 3D polarization vision chart
Apply to binocular balance test and stereo visual function test
- Optotypes can be changed randomly
- Result-speaking function
When the patient answers with direction keys on the controller, the system will show the result by voice.



CE FDA

Both table stand and wall mount are as standard accessories

SPECIFICATIONS

LED monitor.....23-inch TFT Panel, 1920x1080 pixels, 250CD/m²
Optotype.....E, C(Landolt),Alphabet, Numbers, Children
Visual unit.....Decimal,Snellen(ft,20/20),Snellen(m,6/6)
Dis.....1M~6M (in steps of 0.1m)
Mask type.....Single mask, Horizontal mask, Vertical mask,
Red/Green mask, White/Black mask,
Contrast adjustment, Red/Green adjustment

Standard accessories



Polarized filters



Remote



Wall mount




U disc



Power Adapter

Special test masks.....Ishihara Test, Cross cylinder,
Astigmatism Test, Cross Grid,
Fixation point Test, Red/Green Panel,
Worth Test, Red/Green Balance,
Traffic light, Coincidence Test, Stereo Test,
Schober Test, ETDRS
Other functions.....Ramdomization Optotypes Playing function,
R/G setting function, Video playing function,
Screen-saver function, standby function



ZN-185
ZN-215

LED VISUAL CHART
LED VISUAL CHART

18.5"
21.5"HD



ZN-185 CE FDA



ZN-215 CE FDA

Both table stand and wall mount are as standard accessories

SPECIFICATIONS

LED monitor.....ZN-185:18.5-inch TFT Panel, 1366x768 pixels
ZN-215:21.5-inch TFT Panel, 1920x1080 pixels
Brightness.....250CD/m²
Optotype.....E, C(Landolt),Alphabet, Numbers, Children
Visual unit.....Decimal,Snellen(ft,20/20),Snellen(m,6/6)
Distance.....1M~6M (in steps of 0.1m)
Mask type.....Single mask, Horizontal mask, Vertical mask,
Red/Green mask, White/Black mask,
Contrast adjustment, Red/Green adjustment
Special test masks.....Ishihara, Cross cylinder, Astigmatism,
Cross Grid, Fixation point Test, Red/Green Panel,
Worth Test, Red/Green Balance, Traffic light,
Coincidence Test, Stereo Test, Schober Test, ETDRS
Other functions.....Ramdomization Optotypes Playing function,
R/G setting function, Video playing function,
Screen-saver function, standby function

Standard accessories



Wall mount



Power adapter



Remote

