



Miis

Medical Imaging Innovation Solution Provider

DIGITAL HAND-HELD

DIAGNOSTIC SET

Easy take , Easy Record,
Easy Communication



horus
SCOPE




- Multi-Function Interchange Lens
Portable for mobile medicine
- Full HD digital image record
- Real Time Dual-Display on PC and Monitor
Wireless transmit data to the mobile device.



PHYSICIAN



DIAGNOSTIC
STATION



SPECIALIST
NETWORK



TELEMEDICINE
& EMR



MOBILE
MEDICINE



REMOTE
MEDICINE



Wi-Fi
compatible


- AV cable
compatible

fili 3 hours
operation


 USB
compatible

DIGITAL HAND-HELD
DIAGNOSTIC SET


**Fundus
Scope 25°/40°**



Otoscope





Dermoscope




**Digital
Endoscope**

CT



**General
Scope**



Specification

Eye-fundus camera	Non-mydratic, 40 degree FOV, +/-20D
Otoscope	Wide focus range 5mm-50mm
Dermoscope	10mm, Polarization/Contact modes
General lens	f=4mm, F/1.5
Resolution	Full HD True 1920 X 1080
Light Source	Natural white LED
LCD Display	3.5" full color LCD
File Format	JPEG (still picture), H.264 (Video)
Interface	Mini USB and AV output
Memory	Micro SD card (2GB to 32GB)
Power	Li-ion rechargeable battery
External Power	Source: 100-240 VAC, 50/60 Hz
Operating Time	3 Hours
Weight	350 Grams (Approx.)