



neo cosmos



## ■ Leading company in the **global market**



## The Business Philosophy

### ● Human Dignity

Cultivate men of good ability

### ● Customer Focus

Customer-centered strategy

### ● Challenge spirit

A lasting challenge toward the new market

### ● Globalization

The successful realization for the global management

### ● R&D


Research and development should take precedence over the others



## ■ History of Neo Visionce is to Strive For Excellence

- 1993. 07. To establish Poong Gwang Optical Research Laboratory by the assistance of high polymers lesson, technology of Jeon Nam University.
- 1999. 11. Since beginning to export to Singapore, we could achieve a great work of the diversification of export markets.
- 2001. 09. To move the office and change the name of the company to Neo Vision.
- 2002. 10. To be awarded the certificate of SFDA for sales in China and to establish an individual cooperation.
- 2003. 07. EN ISO 9001 certified. (UNK20559)
- 2004. 07. To secure CE certificate. (EN ISO 13488:2000, EN46002:1996 etc.)
- 2005. 05. To apply for a patent of cosmetic color lens material (License No.10-2005-004792) and coating technology.
- 2006. 06. To start processing on commission with Chinese factory.
- 2006. 11. To obtain a patent for Sandwith process. (License No. 10-0647133)
- 2007. 02. KGMP certified by KFDA. (No. MGK-25)
- 2008. 04. To establish the second factory with cast mold system in Wonju city.
- 2008. 07. KFDA approval MPS Neo Plus.
- 2009. 03. Nominated as "2009 The Korea First Class Company" -Herald economy  
Nominated as "2009 BEST INNOVATION Company& Brand"-Sports Seoul.
- 2010. 12. Chairman awarded by Innobiz Association of Korean SMBAs.
- 2011. 06. Awarded as "Promised Exporting Corporation" by Korea SMBA.
- 2012. 05. Begin construction fully automated 3rd factory in Osong city.
- 2013. 05. To establish the third factory with cast mold system in Osong city.
- 2014. 10. To obtain the FDA approval with Color contact lens. (polymacon)
- 2016. 02. To establish Chinese branch in Shanghai.

## Sandwich Process



Color Between Two Layers




### Patented sandwich process technology provides comfort and natural-looking effect for all gorgeous eye color

NEO COSMO is manufactured on the sandwich engineering method, sandwich process is to infiltrate the cosmetic color between the inside and outside layer of the contact lens, which makes Comfortable feeling and prevents decolorization of the lens. Neo Vision obtained the patent with this technology in 2006. (License No.10-0647133)



## Aqua Capsule



Cross sectional view of lens (Thickness: 0.10mm)

### Aqua Capsule technology for all day long moist

Neo Vision's own technology – Aqua Capsule is the new material combined with water factors of high-bioaffinity which forms moisture layer on the lens surface. This superior technology brings all day long moist and best comfort with its water contents on lens surface.

## International Certificates







# NEO ONE DAY


WE ALWAYS RESPECT AND CHERISH THE WELLNESS AND BEAUTY OF YOU.  
TO SATISFY YOUR HIGH LEVELED DEMAND, WE RELEASED NEW PRODUCTS.  
ENJOY NEO VISION'S BRAND NEW ONE DAY USE SOFT CONTACT LENS.

## neo 1 DAY cosmo



Give a special gift to yourself. Start your days with fresh feeling.  
**NEO COSMO** will awake your potential charm inside of you.  
Joyful changes are coming every day with Everyday.


B.C(mm)	8.6
Dia(mm)	14.2
Power Range	0.00 ~ -6.00 (0.25 Step) / -6.50 ~ -10.00 (0.50step)
Materials	METHAFILCON A
Water Content	55%
Color	Blue, Green, Violet, Brown, Gray, Aqua, Honey, Black
Manufacturing Process	Cast Molding



BROWN/ D014




BROWN/ D024



BLACK/ D029



BLUE/ D031




GREEN/ D032



BROWN/ D034



GRAY/ D035



AQUA/ D036



HONEY/ D038

## neo 1 DAY clean



Any time, any place, widen your sight with easy use.  
**Neo Clean's** UV blocker will protect your eyes against intense sunlight.  
**Neo Clean** is the best friend of your clear and healthy eyes.  
Feel the freshness every day.

B.C(mm)	8.6
Dia(mm)	14.2
Power Range	0.00 ~ -6.00 (0.25 Step) / -6.50 ~ -10.00 (0.50step)
Materials	METHAFILCON A
Water Content	55%
Color	Blue Tinted
Manufacturing Process	Cast Molding



30 pcs



# Style A

## N22

2 Color System



BLUE / N221



GREEN / N222



VIOLET / N223



BROWN / N224



GRAY / N225



AQUA / N226



HAZEL / N227



HONEY / N228

## N33

3 Color System



BLUE / N331



GREEN / N332



VIOLET / N333



BROWN / N334



GRAY / N335



AQUA / N336



HAZEL / N337



HONEY / N338

## N42

4 Color System



BLUE / N421



GREEN / N422



BROWN / N424



GRAY / N425

# neo cosmo style A Amore



### SPECIFICATION


B.C.(mm)	8.6
Dia(mm)	14.2
Power Range	0.00 ~ -6.00 (0.25 Step) / -6.50 ~ -10.00 (0.50step)
Materials	2-HEMA
Water Content	38% / 45%
Color	Blue, Green, Violet, Brown, Gray, Aqua, Hazel, Honey
Product Life	Monthly / Quarterly / Yearly
Manufacturing Process	Semi Cast Molding

A deep and rich color shines Amore through in your eyes.  
It expresses younger's sense  
and provide them with the feeling of love.



# NEO PLUS

## MULTI-PURPOSE CONTACT LENS SOLUTION



multi  
purpose **SOLUTION**

NEO-PLUS, multi-purpose contact lens solution has the functions such as; disinfection, cleaning & rinsing, protein removal and preservation which don't have any sore eyes when in use.

capacity ■ 60ml ■ 130ml ■ 360ml




# SOLUTION

## NEO PLUS


NEO-PLUS, multi-purpose contact lens solution has the functions such as; disinfection, cleaning & rinsing, protein removal and preservation which don't have any sore eyes when in use.



**STEP 1**  
Disinfection



**STEP 2**  
Cleaning & Rinsing



**STEP 3**  
Protein Removal



**STEP 4**  
Preservation



•60ml

•130ml


•360ml

## SILICOS

Silicone Hydrogel Lens

*'New evolution of color contact lens'*

Neo Vision's color silicone hydrogel lens is made of elastic material which has higher oxygen permeability compared to conventional hydrogel lenses and has excellent comfort for long time use.



Aqua Capsule

55% Water

Cross sectional view of the lens  
(Thickness 0.10mm)

B.C(mm)	8.6
Dia(mm)	14.2
Power Range	0.00 ~ -6.00 (0.25 Step) / -6.50 ~ -10.00 (0.50step)
Materials	Silicone Hydrogel
Water Content	55%
Color	Blue, Green, Brown, Gray, Aqua, Honey, Black
Manufacturing Process	Cast Molding



**EUROTECH OPTICAL**  
TEL: 012 370 4175 / 1339 / 4690 / 3383 / 3951 / 3470  
email : [sales@eutotechoptical.com](mailto:sales@eutotechoptical.com) [www.eurotechoptical.com](http://www.eurotechoptical.com)  
For the best prices & Quality on Optical Frames & Contact Lenses

**EUROTECH OPTICAL**  
TEL: 012 370 4175 / 1339 / 4690 / 3383 / 3951 / 3470  
email : [sales@eutotechoptical.com](mailto:sales@eutotechoptical.com) [www.eurotechoptical.com](http://www.eurotechoptical.com)  
For the best prices & Quality on Optical Frames & Contact Lenses