

# L4000 Slit Lamp Microscope

## NIPPON SLIT LAMP

### SPECIFICATIONS:

- Microscope: **Stereoscopic**
- Magnification: 10x, 16x, 25x
- Eyepiece: 12.5x
- Slit Width: **0-12mm**
- Slit Height: **1-12mm**
- Slit Apertures:  $\phi 12$ 、 $\phi 9$ 、 $\phi 5$ 、 $\phi 3$ 、 $\phi 1$ 、 $\phi 0.2$  (mm)
- Filters: **Heat Absorption, Grey, Cobalt Blue, Redfree**
- Slit Angle: **0°-180°**
- Slit Illumination: **12V 30W, Halogen lamp**
- Power source: **A.C. 110V or 220V, 50/60Hz**
- Net dimension: **660×285×325mm**
- Net weight: **12Kg**
- Package dimensions: **690×530×235mm**
- Gross weight: **15kg**

