## assistLine

Versatile and space saving.



**OPHTHALMIC SURGERY** 



**AESTHETIC SURGERY** 



DENTAL **SURGERY** 



**FNT** SURGERY



MICRO SURGERY



ORAL & MAXILLOFACIAL

Order No. 24000015

#### **ADVANTAGES**

- Rugged safety chassis with 5-arm aluminum pedestal
- Double casters for extremely smooth movement
- Height adjustment via manually triggered gas pressure spring
- Flexible height adjustment from 54.5 to 73.5 centimeters
- Rotating and tilting anatomically correct saddle seat, comfortable and individually adjustable
- Maximum comfort level for your assistants, to prevent pain and fatigue while supporting the blood circulation
- High-quality, comfortable and easy-care upholstery to meet even the most stringent demands and hygiene requirements
- Seat upholstery in standard color anthracite (Color chart Tundra skai)
- Upholstery optionally available in a wide variety of colors
- Supports maximum loads up to 130 kilograms



Manual gas pressure spring adjustment.

#### **SPECIFICATIONS**

Lifting Column

Seat

- Hand operated

190 mm - Seat adjusting

- Height adjustment 54.5-73.5 cm

Chassis - Five legged chassis ø 585 mm

Upholstery

- Wheels

7,5 kg

Weight Saddle seat

+2°/+16°

Tundra skai

ø 65 mm Standard color anthracite

35 x 42.5 x 18 cm Subject to technical modifications.







### **assistTrend**

**ACCESSORIES** 

### assistPro F

**ACCESSORIES** 

# assistPro H, easyPro, assistLine ACCESSORIES



ARM REST PAIR incl. traverse
Order No. 33000000



DOUBLE RIMMED CASTOR 50mm, elektrically conductive Order No. 33510039



ARM REST PAIR incl. traverse
Order No. 33510038



ELECTRICALLY CONDUCTIVE UPHOLSTERY In color black Order No. 30370223



DOUBLE RIMMED CASTOR 50mm, elektrically conductive Order No. 33510039



DOUBLE RIMMED CASTOR 65mm, with braking function Order No. 33510040



DOUBLE RIMMED CASTOR 65mm, with braking function Order No. 33510040



DOUBLE RIMMED CASTOR 50mm, elektrically conductive Order No. 33510039



DOUBLE RIMMED CASTOR 65mm, with braking function Order No. 33510040

ELECTRICALLY CONDUCTIVE UPHOLSTERY In color black Order No. 30370223