

## EQUIPMENT



Flexible, adaptable and sturdy - the backrest moves synchronously with the footrest (55°) as the chair is electronically reclined to a recumbent position. The chair can also be adjusted to the stable "Trendelenburg" Position.



Wheel brakes on the left and right 360-degree rotating wheels hold the chair safely in place.



Flexible, height adjustable neck support pillow as standard version for enhanced patient relaxation.



Ergonomically formed, durable, non-sliding rubber grips.

Synchronized back-/footrest and seat adjustment via hand control.



High capacity, eco-friendly, independent battery box with emergency stop button.



Comfortable upholstery, anthracite F6461197 standard shade.

## SPECIFICATIONS

### Chassis

- Chassis length 890 mm
- Chassis width 528 mm
- Width incl. armrest 650 mm
- Wheel diameter  $\varnothing$  300/125 mm
- Emergency stop button

### Lifting Column

- 3 parts support frame
- Total length 1,710 mm
- Width 520 mm
- Upholstery Height 40/60 mm
- Back/foot adjustment 55°
- Armrests 0°/+180°
- Seatrest adjustment to Trendelenburg position 25°
- Electronic control system per keypad

- Synchronized back-/footrest and seat adjustment

### Battery unit

- Battery box with 2 batteries (lead-gel), pollution free, maintenance-free, cadmium-free

2 x 12 V 7.2 Ah

- Universal battery charger:

Input: 100-240 V, Output: 24 V

Max.Load 135 kg

Weight 64 kg

Upholstery colors Tundra Skai

Standard color anthracite F6461197

Subject to technical modifications.

