

## EQUIPMENT



Central brake acting on 4 wheels.



Both seat height adjustment and armrest height adjustment are effected electronically via the foot control element.



Ergonomically shaped push handles for easy maneuvering.



The individually adjustable armrests can be adjusted under sterile conditions via the foot keypad and the one-hand release.



"Dynamic"® seating system: our unique dynamic seating system supports correct sitting posture while protecting your back.



Optional: surgiForce with Combi Seat. **No. 26000001**

Individual stool adjustment.

The optionally adaptable foot switch plate ensures the optimized layout of your foot control elements even during alternating uses, thus obviating the need for time-consuming repositioning maneuvers. **No. 33170000**



Exemplary application with foot switch.

## SPECIFICATIONS

### Lifting Column

- Lifting height of the seat 200 mm
- Lifting height of the armrests 135 mm
- Electromotive adjustment via the foot control elements

### Chassis

- w x l x h max. (incl. armrest) 742 x 861 x 943 mm
- Wheels ø 100 mm
- 5. Wheel ø 60 mm
- Central brake acting on all four wheels
- Emergency stop button

### Battery unit

- Lead-gel battery, 2 ea. 12 V, 2.9 A
- Universal battery charger:
- Input: 100-240 V
- Output: 24 V

### Seat/backrest

- Min. sitting height 640 mm
- Max. sitting height 840 mm
- Seat/Back adjusting -6°/+3°

### Weight

74 kg

### Max Load

150 kg

### Upholstery colors

Stamskin Tpp

### Standard color

black

Subject to technical modifications.

