


EQUIPMENT



"Dynamic"® seating system: the dynamically synchronized lumbar-support seating mechanism can be individually adjusted to every user, thus promoting healthy seating behavior. During frequent switching between sitting and standing work positions, every change of position is actively supported.



Optional: one-hand height release for direct and high-precision all-around adjustment.

No. 33100000



Central braking system effective on 4 wheels.



Foot keypad panel-operated electro-motive drives for height adjustment.

Optional special accessories: Storage module incl. cable guide for microscope/phaco foot switches.

No. 32070000



*Optional special accessories: Adaptable foot-switch plate for your control elements on 5 casters. To adjust the foot switch together with the seat to best suit your individual position. Useful, convenient, and time-saving. **No. 31050003***

SPECIFICATIONS

Lifting Column

- Drive 3,000 N IPX4
- Lifting height 200 mm

Chassis

- w x l 470 x 535 mm
- Wheels ø 100 mm
- 5. Wheel ø 65 mm
- Central brake acting on all four wheels
- Emergency stop button

Weight 45 kg

Max Load 150 kg

Battery unit

- Battery box with 2 batteries

(lead-gel), pollution free, maintenance-free, cadmium-free
2 x 12 V 7.2 Ah

- Universal battery charger:

- Input: 100-240 V
- Output: 24 V

Seat/backrest

- Min. sitting height 553 mm
- Max. sitting height 753 mm
- Seat/Back adjusting -6°/+3°

Upholstery colors Stamskin top

Standard color black

Subject to technical modifications.

