

SPECIFICATIONS

Description	Medical Ultrasonic Cleaner	
Model	CS8.0	
Tank Capacity	8.0L / 2.1 gallon (US)	Max. 7.0 L / 1.8 gallon (US)
		Min. 3.6 L / 0.95 gallon (US)
Tank Size	33.1x23.8x11.6 cm / 13.0"x9.4"x4.6"	
Longest Item Fits inside Tank	37.0 cm / 14.4"	
Power Supply	450W (AC 100~120V 50/60Hz)	Ultrasound Power: (68-75)Wx3=(204-225)W
	450W (AC 220~240V 50/60Hz)	3 ceramic heaters: (70-75)Wx3=(210-225)W
Timer settings	☞ 5min → 10min → 15min → 20min → 25min → 30min ☞	
Temperature settings	☞ 40℃ → 45℃ → 50℃ → 55℃ → 60℃ ☞	
	☞ 104°F → 113°F → 122°F → 131°F → 140°F ☞	
Drainage	Drainage valve	
Ultrasonic Frequency	46KHz ± 3KHz	
Tank Material	Stainless Steel SUS304	
Housing Material	Aluminum+ABS	
Net Weight	14.3 kg / 31.5lb	
Gross Weight	17.1 kg / 37.7lb	
Unit Size	46.6x35x22.9 cm / 18.3"x13.8"x9.0"	
Inner Carton Size	52.5x39x27.5 cm / 20.7"x15.3"x10.8"	
Qty per Master Carton	1 pc / ctn	
Master Carton Size	54x40.5x30.5 cm / 21.3"x16.0"x12.0"	