

Neo Cosmo Solution Neo Plus



Multi-purpose contact lens solution

Rating: Not Rated Yet
[Ask a question about this product](#)

Description

Multi-purpose contact lens solution

[Click here to view the pdf catalogue](#)

NEO-PLUS, multi-purpose contact lens solution has the functions such as; disinfection, cleaning & rinsing, protein removal and preservation which don't have any sore eyes when in use.

Sandwich Process

Color Between Two Layers

Safe, Secure and Sandwiched!

Beauty zone

Optical zone

External polymer

Internal coating

Patented sandwich process technology provides comfort and natural-looking effect for all gorgeous eye color

NEO COSMO is manufactured on the sandwich engineering method, sandwich process is to infiltrate the cosmetic color between the inside and outside layer of the contact lens, which makes Comfortable feeling and prevents decolorization of the lens.

Neo Vision obtained the patent with this technology in 2006.
(License No.10-0647133)

[SANDWICH] TRADEMARK [SANDWICH] PATENT

