

Shinva Tunnel Washer-Disinfector DC Series



Washer-Disinfector

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

[Click here to view the pdf catalogue](#)

Introduction

Tunnel washer-disinfector is internationally-advanced product which is researched aim at middle and large-scale hospitals. It has advantages of reliable quality, high work load and high efficiency and can be widely used for hospital reusable articles such as operating instrument, minimally invasive instrument, corrugated and anesthesia hose, glass ware, dental handpiece, baby product and general medical instrument. It is also applicable for hospital CSSD.

Modular design, compact structure, reliable performance, convenient installation It adopts modular design and provides multiple combinations to meet customer requirements. The wide of the washer-disinfector is only 1200mm which provides convenient installation and furthest reduces the expenses and time of installation.

Product features

- **More efficient**
Advanced design concept which makes the running speed 30% faster than traditional products.

- **Lager capacity**

The lager capacity which can afford 15 DIN trays one time gives it 50% more than traditional products.

- **Better washing result**

Well-designed spray arm with flexible rotary system and perfect connection between washing rack and chamber provide a better washing result than traditional products.

- **Easy-to-use**

To operate this washer-disinfector, simply touch the 10.6 inches touch-screen.

- **Energy saving**

New-designed pipeline system and dead-zone free structure save 25% water.

- **Reliable and sanitary inner-chamber delivery system**

The chamber equipped with stainless steel wheels which are dead-zone free. All transfer executive components are installed out of chamber to protect the instrument from secondary contamination.

- **Better compatibility**

The washing racks can also work with Shinva Rapid-A-520 or Super6000 washer-disinfector.

- **Humanized window design and good maintainability**

The chamber adopts glass-faced sealing door and large-sized window on side of it. For maintenance purpose, all pipelines are fixed on one side of the washer-disinfector with open type door.

- **High-efficient hot-wind drying system**

Two-stage drying system with sieve-style air inlet dramatically improves drying result and efficiency. It also solves the drying problems caused by temperature changes.

- **Diversified carrying system**

For your convenience, we offer manual carry-trolley and auto washing rack carrying system to increase efficiency and reduce labor intensity.